



Promotional article

Maxi-Smart

For over twenty years, Cattani has been doing research and is currently up-to-date with the latest technology.

The primary purpose of research is to constantly improve our products: building machines that are increasingly flexible, safer, less expensive to buy and, most importantly, to run and service.

In addition, Cattani has followed a personal evolution path, based on scientific and technological research, with focus on sustainable growth and on solutions that allow us to save on raw materials and energy.

Below we are presenting two new products: Maxi-Smart with box e Maxi-Smart without box.

Maxi- Smart is a dental surgery suction system, recommended for 15 treatment units running simultaneously.

Description:

- the centrifugal separator replaces both the separator tank and draining pump, thus reducing dimensions and an important part of periodic maintenance;
- the amalgam separator is assembled with the centrifugal separator on top of the suction system and next to the secretion filter, making these two parts very easy to reach for routine maintenance;
- the suction unit is covered and soundproofed, residual vibration is 68 dB (A) (in the version without box the sound pressure is 77dB(A));
- the control unit contains both the inverter and a small computer;
- the programme includes optional head adjustment,
- automatic flow rate adjustment according to use,
- and remote danger and control warnings via both the Web and App;
- Maxi-Smart is the first large system produced as a single unit;





- Maxi-Smart is absolutely unique at global level: it is the first large surgery suction system that can control flow rate as well as head;
- with Maxi-Smart simultaneous use no longer means reduced performance: as demand increases, flow rate automatically, instantly increases.

Augusto Cattani

